

# 730050 - ColorCheck Color Vision Recording Form

## Number Color Plates

| Plate Number              | Normal | Mild Defect        | Severe Defect   | Patient's Response |
|---------------------------|--------|--------------------|-----------------|--------------------|
| 1 - Demonstration Plate   | 35     | 35                 | 35              |                    |
| 2                         | 9      | Nothing Seen       | Nothing Seen    |                    |
| 3                         | 2      | Nothing Seen       | Nothing Seen    |                    |
| 4                         | 73     | 28 or Nothing Seen | Nothing Seen    |                    |
| 5                         | 34     | 81 or Nothing Seen | Nothing Seen    |                    |
| 6                         | 7      | Nothing Seen       | Nothing Seen    |                    |
| 7                         | 92     | Nothing Seen       | Nothing Seen    |                    |
| 8                         | 35     | Nothing Seen       | Nothing Seen    |                    |
| 9                         | 57     | Nothing Seen       | Nothing Seen    |                    |
| 10                        | 3      | Nothing Seen       | Nothing Seen    |                    |
| 11                        | 74     | 21 or Nothing Seen | Nothing Seen    |                    |
| 12                        | 95     | Nothing Seen       | Nothing Seen    |                    |
| 13                        | 5      | 3 or Nothing Seen  | Nothing Seen    |                    |
| 14                        | 4      | Nothing Seen       | Nothing Seen    |                    |
| 15                        | 93     | 98 or Nothing Seen | Nothing Seen    |                    |
| 16                        | 75     | (7)5 Mild Defect   | 5 Severe Defect |                    |
| 17 - Tritan (Blue/Yellow) | 3      | Nothing Seen       | Nothing Seen    |                    |
| 18 - Tritan (Blue/Yellow) | 57     | Nothing Seen       | Nothing Seen    |                    |

# 730050 - ColorCheck Color Vision Recording Form

## LEA SYMBOLS® Color Plates

| Plate Number             | Normal  | Color Vision Defect | Patient's Response |
|--------------------------|---------|---------------------|--------------------|
| 19 - Demonstration Plate | ♥ □ ◡ ○ | ♥ □ ◡ ○             |                    |
| 20                       | □ ○     | □                   |                    |
| 21                       | ♥ □     | ♥                   |                    |
| 22                       | □ ♥     | □                   |                    |
| 23                       | ○ □     | ○                   |                    |
| 24                       | ♥ □     | Nothing Seen        |                    |
| 25                       | □ ◡     | Nothing Seen        |                    |
| 26                       | ○ □     | Nothing Seen        |                    |