

CHART PROJECTOR

CP-500



Rexxam

CHART PROJECTOR

CP-500



50 charts and various masks

Landolt, Letter, Child, Snellen, 50 kinds of charts, and Horizontal line, Vertical line, one chart, various masks can be selected by the remote control box.

Fast, silent chart disk rotation

It moves to the following chart in a maximum of 0.3 seconds.

2 kinds of programs can be selected. Each program can store 48 kinds (maximum) of chart and mask.

Various Switches

- 1 Lamp switch
- No use
- 4 Chart/Mask switch One chart switch
- 6 Horizontal line switch
- 8 Back switch
- 3 Program switch (1,2)
- 6 Vertical line switch
- 9 Next switch

SPECIFICATIONS	
Number of chrats and masks	50 charts and 38 masks
Luminance	min.80-max.320Cd/m2(200Cd/m2 at 5m)
Chart and mask selection speed	0.3seconds
Program	2 kinds(48 charts or masks)
Tilting angle	±10°
Auto power off	After 5 min.
Lamp	White LED
Projection distance	2~7 m
Projection magnification	30 x (at 5m)
Projection size	335mm x 222mm (at 5m)
Dimentions	W248xL336xH234mm
Weight	4.5kg
Power source	100~240V
Power consumption	30VA
STANDARD ACCESSORIES	
REMOTE CONTROL BOX, SCREEN, TEST SCALE	E SHEET, DUST COVER, SPARE FUSE

OPTIONAL ACCESSORIES

TABLE BASE STAND, WALL MOUNT, CHART PROJECTOR ARM, POLA FILTER GLASS, R&G GLASS, DEVIATOR

The brightness of the projector can be changeable by the control switch according to the brightness of the operation room.



A type chart (B type chart is also available.



5 2 8 0.1











8 4 9 5 0 0.9 0 3 7 2 5 1.0

35792 0.8

703680.8



5 2 8 3 9 0.9

496571.0















Rexxam Co., Ltd.

Contact

Eye-care Instruments Sales Dept. Tokyo Office

2-4-2 Kandatsukasa-machi, Chiyoda-ku, Tokyo, 101-0048, Japan TEL:81-3-6262-9471 FAX:81-3-6262-9472

E-mail:eye@rexxam.co.jp

URL:http://www.rexxam.co.jp URL:http://www.shin-nippon.co.jp

Manufacturer



TOWA MEDICAL INSTRUMENTS Co., Ltd.

1698-2 Oaza Nakano Nakano-shi, Nagano-ken, 383-0013 Japan