

GAME CHANGER IN

CONTRAST SENSITIVITY TESTING

CONTRAST SENSITIVITY **CAMBLOBS** TEST



Contrast Sensitivity testing is a vital part of testing functional vision. Reductions in contrast sensitivity can indicate the onset or progression of many diseases and severely complicate activities of daily living.

Traditionally, testing contrast sensitivity has been difficult and very time consuming to perform. In some instances, test/retest repeatability and interpretation and communication of results are very discouraging.

CAMBLOBS™ is a new test that changes everything about gathering precise contrast sensitivity testing results. Each test is a single use test where patients test themselves by marking an “x” on each contrast target they see. Testing is intuitive; it only takes 2-3 minutes to complete the test and scoring the results is very easy and takes only a minute. The test sheet itself then becomes the record for the patient file.

Contact Precision Vision for more details about CAMBLOBS and any other complete lines of contrast tests.

A CamBlobs

A CamBlobs2

LogCS	Name/ID:	Date:	L / R / B	LogCS:	3 No.
0.85					
0.9					
0.95					
1.0					
1.05					
1.1					
1.15					
1.2					
1.25					
1.3					
1.35					
1.4					
1.45					
1.5					
1.55					
1.6					
1.65					
1.7					
1.75					
1.8					
1.85					
1.9					
1.95					
2.0					
2.05					

SN: 99999999 © JGRobson 2015

Precision Vision